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1. Chemistry, Manufacturing, and Control (CMC) Information for Human Gene Therapy Investigational New Drug Applications (INDs), Guidance for Industry. US Food and Drug Administration (CDER), January 2020.
2. Kato, Y.; Sasaki, M.; Hashimoto, T. Separation of DNA Restriction Fragments by High-Performance Ion Exchange Chromatography. *J. Chromatogr.* 1983, 265, 342–346.
3. Kato, Y.; Yamasaki, Y.; Onaka, A.; Kitamura, T.; Hashimoto, T. Separation of DNA Restriction Fragments by High-Performance Ion Exchange Chromatography on a Non-Porous Ion Exchanger. *J. Chromatogr.* 1989, 478, 264–268.
4. Muller, W. Fractionation of DNA Restriction Fragments with Ion-Exchangers for High-Performance Liquid Chromatography. *Eur. J. Biochem.* 1986, 155, 203–212.
5. Shapiro, J. T.; Stannard, B. S.; Felsenfeld, G. Binding of Small Cations to Deoxyribonucleic Acid. Nucleotide Specificity. *Biochemistry* 1969, 8(8), 3233–3241.
6. Melchior, W. B. Jr.; Von Hippel, P. H. Alteration of the Relative Stability of dA-dT and dG-dC Base Pairs in DNA. *Proc. Nat. Acad. Sci. USA.* 1973, 70(2), 298–302.
7. Yang, H.; Koza, S.; Chen, W. Anion-Exchange Chromatography for Determining Empty and Full Capsid Contents in Adeno-Associated Virus. Waters Application Note, [720006825 < https://www.waters.com/nextgen/us/en/library/application-notes/2020/anion-exchange-chromatography-for-determining-empty-and-full-capsid-contents-in-adeno-associated-virus.html >](https://www.waters.com/nextgen/us/en/library/application-notes/2020/anion-exchange-chromatography-for-determining-empty-and-full-capsid-contents-in-adeno-associated-virus.html) , 2021.
8. Yang, H.; Koza, S.; Chen, W. Plasmid Isoform Separation and Quantification by Anion-Exchange Chromatography (AEX). Waters Application Note, [720007207 < https://www.waters.com/nextgen/us/en/library/application-notes/2021/plasmid-isoform-separation-and-quantification-by-anion-exchange-chromatography-aex.html >](https://www.waters.com/nextgen/us/en/library/application-notes/2021/plasmid-isoform-separation-and-quantification-by-anion-exchange-chromatography-aex.html) , 2021.

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© 2021 Waters Corporation. All Rights Reserved.