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1. Inoue, K.; Obara, R.; Hino, T.; Oka, H. Development and Application of an HILIC-MS/MS Method for the Quantitation of Nucleotides in Infant Formula. *J. Agricultural and Food Chemistry* 2010, 58(18), 9918–9924.
2. Logotheti, M.; Theochari, K.; Kostakis, M.; Pasiadis, I.; Thomaidis, N. Development and Validation of a HILIC–UV Method for the Determination of Nucleotides in Fish Samples. *Food Chemistry* 2018, 248, 70–77.
3. Mateos-Vivas, M.; Rodriguez-Gonzalo, E.; Domingues-Alvarez, J.; Garcia-Gomez, D.; Carabias-Martinez, R. Determination of Nucleosides and Nucleotides in Baby Food by Hydrophilic Interaction Chromatography Coupled to Tandem Mass Spectrometry in the Presence of Hydrophilic Ion-Pairing Reagents. *Food Chemistry* 2016, 211, 827-835.
4. Zhao, H.; Wang, X.; Li, H.; Yang, B.; Yang, H.; Huang, L. Characterization of Nucleosides and Nucleobases in Natural Cordyceps by HILIC-ESI/TOF/MS and HILIC-ESI/MS. *Molecules* 2013, 18(8), 9755–9769.
5. Chen, P.; Li, W.; Li, Q.; Wang, Y.; Li, Z.; Ni, Y.; Koike, K. Identification and Quantification of Nucleosides and Nucleobases in Geosaurus and Leech by Hydrophilic-Interaction Chromatography. *Talanta* 2011, 85(3), 1634–1641.
6. Jung, M.; Lauber, M. MaxPeak高性能表面(HPS)技术在改善灵敏度和动态范围方面的表现：关于核苷酸检测的案例研究. 沃特世应用纪要, [720007053ZH](https://www.waters.com/nextgen/us/en/library/application-notes/2020/demonstrating-improved-sensitivity-and-dynamic-range-with-maxpeak-high-performance-surfaces-hps-technology-a-case-study-on-the-detection-of-nucleotides.html) <<https://www.waters.com/nextgen/us/en/library/application-notes/2020/demonstrating-improved-sensitivity-and-dynamic-range-with-maxpeak-high-performance-surfaces-hps-technology-a-case-study-on-the-detection-of-nucleotides.html>>, 2020.

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- [ACQUITY Premier系统](https://www.waters.com/waters/nav.htm?cid=135077739) <<https://www.waters.com/waters/nav.htm?cid=135077739>>
- [Xevo TQD三重四极杆质谱仪](https://www.waters.com/134608730) <<https://www.waters.com/134608730>>

- [MassLynx MS软件 <https://www.waters.com/513662>](https://www.waters.com/513662)
- [Atlantis Premier色谱柱 <https://www.waters.com/waters/nav.htm?cid=513211>](https://www.waters.com/waters/nav.htm?cid=513211)

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© 2021 Waters Corporation. All Rights Reserved.