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---

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1. Erickson, BE.(2021).C&En News.99(31)."Delta-8-THC Craze Concerns Chemists - cen.acs.org."30 Aug. 2021, <https://cen.acs.org/biological-chemistry/natural-products/Delta-8-THC-craze-concerns-chemists/99/i31> < <https://cen.acs.org/biological-chemistry/natural-products/Delta-8-THC-craze-concerns-chemists/99/i31>> .
2. Hudalla, C. (2021). "We Believe in Unicorns (and Delta-8)" .The Cannabis Scientist. <https://thecannabisscientist.com/testing-processing/we-believe-in-unicorns-and-delta-8> < <https://thecannabisscientist.com/testing-processing/we-believe-in-unicorns-and-delta-8>> .
3. Webster, GRB., Sarna, LP., Mechoulam, R. US Patent Application Publication.Pub. <https://patents.google.com/patent/US20040143126A1/en> < <https://patents.google.com/patent/US20040143126A1/en>> .
4. Golombek, P., *et al.* (2020)."Conversion of Cannabidiol (CBD) into Psychotropic Cannabinoids Including Tetrahydrocannabinol (THC): A Controversy in the Scientific Literature."*Toxics* 8(2).
5. Kiselak, T. D., *et al.* (2020)."Synthetic Route Sourcing of Illicit At Home Cannabidiol (CBD) Isomerization to Psychoactive Cannabinoids Using Ion Mobility-Coupled-LC-MS/MS."*Forensic Sci Int* 308: 110173.
6. Meehan-Atrash, J., and I. Rahman (2022)."Novel  $\Delta$ 8-Tetrahydrocannabinol Vaporizers Contain Unlabeled Adulterants, Unintended Byproducts of Chemical Synthesis, and Heavy Metals."*Chemical Research in Toxicology* 35(1): 73–76.
7. Watanabe, K., *et al.* (2007)."Conversion of Cannabidiol to  $\Delta$ 9-Tetrahydrocannabinol and Related Cannabinoids in artificial gastric juice, and their pharmacological effects in mice."*Forensic Toxicology* 25(1): 16–21.
8. Helander, A., *et al.*(2022)."Analytical and medico-legal problems linked to the presence of delta-8-tetrahydrocannabinol (delta-8-THC): Results from urine drug testing in Sweden."*Drug Testing and Analysis* 14(2): 371–376.
9. Gorenstein, MV., Li, JB., Van Antwerp, J., Chapman, D., Detecting Coeluting Impurities by Spectral Comparison.<https://www.waters.com/webassets/other/SpectralContrastAlgorithm.PDF> <

<https://www.waters.com/webassets/other/SpectralContrastAlgorithm.PDF>> .

10. Lander, N., [Get Empowered: Creating, Searching, and Restoring UV Libraries | Tip #144, Creating a UV Library - Waters](https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-144-creating-a-uv-library) <<https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-144-creating-a-uv-library>> .
11. Lander, N., [Get Empowered: Creating, Searching, and Restoring UV Libraries | Tip #145, Manually Searching a UV Library - Waters](https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-145-manually-searching-a-uv-library) <<https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-145-manually-searching-a-uv-library>> .
12. Lander, N., [Get Empowered: Creating, Searching, and Restoring UV Libraries | Tip #146, Performing an Automated Search in a UV Library - Waters](https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-146-performing-an-automated-search-in-a-uv-library) <<https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-146-performing-an-automated-search-in-a-uv-library>> .
13. Lander, N., [Get Empowered: Creating, Searching, and Restoring UV Libraries | Tip #147, Back up and Restore PDA Libraries - Waters](https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-147-back-up-and-restore-pda-libraries) <<https://blog.waters.com/get-empowered-creating-searching-and-restoring-uv-libraries-tip-147-back-up-and-restore-pda-libraries>> .
14. Lander, N., [Get Empowered: Tip #251 | Peak Purity - Waters](https://blog.waters.com/empower-251-peak-purity?preview_id=11242) <[https://blog.waters.com/empower-251-peak-purity?preview\\_id=11242](https://blog.waters.com/empower-251-peak-purity?preview_id=11242)> .
15. Lander, N., [The Math Behind Peak Purity Calculations | Waters Blog](https://blog.waters.com/empower-252-peak-purity?preview_id=11295) <[https://blog.waters.com/empower-252-peak-purity?preview\\_id=11295](https://blog.waters.com/empower-252-peak-purity?preview_id=11295)> .

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[ACQUITY UPLC PDA 検出器](https://www.waters.com/514225) <<https://www.waters.com/514225>>

[ACQUITY QDa 質量検出器](https://www.waters.com/134761404) <<https://www.waters.com/134761404>>

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© 2024 Waters Corporation. All Rights Reserved.

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